

Attaching straps to the masks: OLSON

Stop at Step 8.

- 8 With the fronts of the fabric facing each other, sew along the ~~entire perimeter of the mask.~~



1. Do not sew the entire perimeter of the mask. Sew perimeter but leave side 4 open on both sides.



2. Invert the mask right side out.
3. Fold raw edges (side 4) on itself so the raw edges are tucked inside the mask. Pin to keep in place
4. Insert in one inch of unfinished edge of one strap into one corner between the front and back layers. Topstitch and backstitch catching the innermost edge of the end between the layers. Repeat on all corners.
5. Stitch straps in place at yellow line at seam allowances 1/2 inch, 1/4 and as close to edge possible. Be sure straps are secure and will not fall out.



Thank you! Construction Complete!